

# Beach chair helmet head rest 30130

Helmet type head rest 30130 for shoulder and upper arm procedures with "Beach Chair" back rest. Unique sliding function with double joint mechanism, very wide adjustments in all directions. Fits Promerix 60525 and Practico 30125 back rests.



2-rod

## Technical Specifications

Max. patient weight	180 kg
Safe working load (SWL)	11 kg

---

Weight	4.7 kg
Product Code	100030130

## Contact us

**Merivaara Corp.**

Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)