

## ENT-surgery head rest 18143

Gas assisted narrow head rest with contour moulding for comfort, vertical up and down manipulation.

Fits all Merivaara operating tables with "Beach Chair" back rest. For Promerix and Practico the fixing adapter 60146 is needed.

### Order no.

110143 Head rest 18143 with adapter 60146 100  
018 143 Head rest 18143 100 060 146 Adapter  
for Promerix/Practico



2-rod

### Technical Specifications

Width	242 mm
Height	245 mm
Thickness	40 mm

Max. patient weight	325 kg
Safe working load (SWL)	14 kg
Trendelenburg	-45 °
Reverse Trendelenburg	+25 °
Weight	3.4 kg
Product Code	110143, 100018143

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)