


Head rest 60140

All purpose gas spring assisted head rest for wide range of surgical applications. Optional mattresses moulded IS or visco elastic VEF. Accessory rail option EUR or US.



Technical Specifications

Width	540 mm
Height	300 mm
Thickness	80 mm
Max. patient weight	325 kg
Safe working load (SWL)	14 kg
Weight	7 kg
Product Code	100060140

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com