


## Head rest 60140

All purpose gas spring assisted head rest for wide range of surgical applications. Optional mattresses moulded IS or visco elastic VEF. Accessory rail option EUR or US.



### Technical Specifications

Width	540 mm
Height	300 mm
Thickness	80 mm
Max. patient weight	325 kg
Safe working load (SWL)	14 kg
Weight	7 kg
Product Code	100060140

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)