


# Head rest adapter 60146

Adapter for Promerix and Practico operating tables for fixing special head rests, e.g. 18150 for ophthalmology, Doro/Mayfield Cranial Positioning System.



## Technical Specifications

Safe working load (SWL)	14 kg
Weight	2.5 kg
Product code	100060146

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)

