


## Push bar 20593

Push bar hand support for maneuvering the operating table and for supporting the end user's hands when activating head rest articulation. Suitable for 60140 and 60142 head rests, Promerix and Practico operating tables. Diameter 16 mm.



### Technical Specifications

Weight	0.8 kg
Product Code	100020593

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)