

Beach chair head rest 30126

Special multi directional head rest for shoulder and upper arm procedures. For Beach Chair back rest for shoulder and upper arm procedures. Unique sliding function with double joint mechanism, very wide adjustments in all directions. Fits to Promerix 60525 and Practico 30125 back rests. Rotates 360 degrees, adjustments 180 degrees. Vertical adjustment: 65 mm

Order no.

145126 Head rest for shoulder arthroscopy 30126 with adapter 60146 100 030 126 Head rest for shoulder arthroscopy





Technical Specifications

--	--

Max. patient weight	225 kg
Safe working load (SWL)	12 kg
Weight	4.6 kg
Product Code	145126, 100030126

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com