

## Head rest 60142

Double articulated gas spring assisted head rest with side rails for Promerix and Practico. Perfect for lateral positioning of the patient. Helps with spinal alignment. Optional mattresses moulded IS or visco elastic VEF. Including accessory rails (EUR/US).



### Technical Specifications

Width	540 mm
Height	300 mm
Thickness	80 mm
Max. patient weight	325 kg
Safe working load (SWL)	14 kg
Weight	9.8 kg
Product Code	100060142

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)