

# Anaesthesia frame with rotation clamp

- stainless steel frame
- stepless height adjustment
- antimicrobial tightening knob
- possible to expand with 1-2 extension arms (AI1002)



## TECHNICAL SPECIFICATIONS

Product code	AI1010 (EU), AI1011 (US), AI1012 (UK)

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)