

# Anaesthesia screen 111

Adjustable height and angle. Rotates 360 degrees. Made of stainless steel. Universal clamp 10300 is needed for fixing. Rod diameter Ø 19 mm

## Order no.

110111	Anaesthesia screen kit, incl. anaesthesia screen 111, 2 x extension arms 116 and 10300 clamp ( for EUR, US and UK rails)
100 000 111	Anaesthesia screen 111
100 000 116	Extension arm for anaesthesia screen 116



## Technical Specifications

Length	590 mm
Width (base)	20 mm
Weight	1 kg
Product Code	110111, 100000111

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)