K-Tek 500 stirrup pair 12812

The K-Tek Comfort Stirrups for patient positioning for a wide range of gynecologic, urologic and laparoscopic procedures. The Stirrups increase surgeons' operating room performance by permitting simple adjustment of lithotomy and abduction intra-operatively. The K-tek Comfort 500 stirrup is designed for patients with a BMI of greater than 29 with a large flexible lateral fin and can support patients up to 227 kg weight capacity.

- Designed to relieve pressure with soft padding that encapsulates the foot.
- K-TEK™ silicone closure straps making cleaning easy.
- Flexible finned boots accomodates various leg sizes.
- Superior access with lithotomy range of motion from -55° to +85° aids in relieving clinican shoulder/arm strain.



Order information

100012812 K-Tek 500

stirrup pair

100012808 K-Tek™

stirrup cart

100012811 K-Tek™

stirrup boot

pad kit

CL1210 Rail clamp,

EUR and UK

100 011 043U Rail clamp,

US

Technical Specifications

Max. patient weight	227 kg
Weight	6.6 kg

Merivaara Corp. 31.10.2025



Movement angle	+-25
Lithotomy Range	+8555 °
Abduction range	+259 °
Product Code	100012812

Contact us

Merivaara Corp. Tarmontie 2-4 15860 HOLLOLA FINLAND

+358 3 3394 611

merivaara@merivaara.com

Merivaara Corp. 31.10.2025