

## Powered leg holders 60482

Powered leg holders 60482 provides excellent patient comfort and precise surgical positioning, the patient can be positioned from the supine position to lithotomy with ease. The VEF-pad reduces pressure. Velcro straps for security. Integrated clamps included. Suitable with Promerix and Smarter Practico operating tables, replacing a leg section. Adjustment ranges: Support arm (vertical) -95 – +70 degrees, supports (horizontal and vertical) 360 degrees.



### Technical Specifications

Length	600 mm
Height	345 mm
Max. patient weight	275 kg
Weight	3.0 kg
Product Code	100060482

## **Contact us**

**Merivaara Corp.**

Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)