

## Thigh holder 12012

For securing the patient in the lateral position, the curvature of the thigh holder ensures patient comfort and stability. Rotating, complete with antistatic pad. Vertical bar 550mm, horizontal bar 184mm. Rod diameter 16 mm. Includes universal clamp 10200 for EUR, UK and US rails.

### Order no.

100 012 012 Thigh holder  
12012 (1pc)  
, includes  
universal  
clamp  
10200 for  
EUR, UK and  
US rails



## Technical Specifications

Length	100 mm
Height	550 mm
Max. patient weight	325 kg
Safe working load (SWL)	40 kg
Weight	1.8 kg
Product Code	100012012

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)