

Thigh holder 12012

For securing the patient in the lateral position, the curvature of the thigh holder ensures patient comfort and stability. Rotating, complete with antistatic pad. Vertical bar 550mm, horizontal bar 184mm. Rod diameter 16 mm. Includes universal clamp 10200 for EUR, UK and US rails.

Order no.

100 012 012 Thigh holder
12012 (1pc)
, includes
universal
clamp
10200 for
EUR, UK and
US rails



Technical Specifications

| | |
|-------------------------|-----------|
| | |
| Length | 100 mm |
| Height | 550 mm |
| Max. patient weight | 325 kg |
| Safe working load (SWL) | 40 kg |
| Weight | 1.8 kg |
| Product Code | 100012012 |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com