

## 90° knee support 138

Knee support with height adjustable, telescopic support post, radiotranslucent (carbon fibre) and padded, attaches to the adapter. Can be fixed on both sides of the table. For orthopedic extension device 19300. Adjustable height 105-475 mm.



### Technical Specifications

Length	500 mm
Height	610 mm
Max. patient weight	225 kg
Safe working load (SWL)	44 kg
Weight	3.6 kg
Product Code	100000138

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)