

Femoral counter traction post 19112

Femoral counter traction post 19112 for orthopaedic extension device 19300. X-ray translucent carbon fibre perineal post. Height adjustable between 250-335 mm.

Including: • Pads • Perineal posts with crank handle



Technical Specifications

Length	270 mm
Height	800 mm
Max. patient weight	225 kg
Safe working load (SWL)	58 kg
Weight	3.1 kg
Product Code	100019112

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com