

Leg rest with padding 19120

Radiotranslucent leg rest, curved support (carbon fibre) with support post and clamp, attaches to the extension bar. Rod diameter Ø 19 mm. Suitable with orthopedic extension device 19300.



Technical Specifications

| | |
|-------------------------|-----------|
| Length | 330 mm |
| Height | 470 mm |
| Max. patient weight | 225 kg |
| Safe working load (SWL) | 29 kg |
| Weight | 1.9 kg |
| Product Code | 100019120 |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com