


## Femoral counter traction post 19112

Femoral counter traction post 19112 for orthopaedic extension device 19300. X-ray translucent carbon fibre perineal post. Height adjustable between 250-335 mm.

Including: • Pads • Perineal posts with crank handle



### Technical Specifications

| Length                  | 270 mm    |
|-------------------------|-----------|
| Height                  | 800 mm    |
| Max. patient weight     | 225 kg    |
| Safe working load (SWL) | 58 kg     |
| Weight                  | 3.1 kg    |
| Product Code            | 100019112 |

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)