

## Short support arm 19313

Short support arm 19313 for supporting healthy leg with stirrup or leg support, e.g. 20118. Clamp 19135 to be ordered separately.



### Technical Specifications

| Length                  | 50 mm     |
|-------------------------|-----------|
| Height                  | 300 mm    |
| Max. patient weight     | 225 kg    |
| Safe working load (SWL) | 58 kg     |
| Weight                  | 1.5 kg    |
| Product Code            | 100019313 |

## Contact us

**Merivaara Corp.**

Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)