


CS Easy Prone Low (chest and pelvis set) CSM-2699-G

- Fully radiolucent
- Fully sealed cover for infection control
- Allows abdomen to hang freely
- Can be flexibly adjusted to suit many different body shapes
- Allows free movement of the operating table
- Unique shape designed to avoid pressure in the armpit, neck and genital area
- Regulatory compliance: CE, FDA, TGA, UKCA
- MD Class: I



Technical Specifications

Chest	
Length	250 mm
Width	500 mm
Height	150 mm
Pelvis	
Length	336 mm
Width	500 mm
Height	150 mm
Weight capacity	150 kg
Product Code	CSM-2699-G
Manufacturer	Care Surgical Limited

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com