

## CS Easy Prone Medium (chest and pelvis set) CSM-2918-G

- Fully radiolucent
- Fully sealed cover for infection control
- Allows abdomen to hang freely
- Can be flexibly adjusted to suit many different body shapes
- Allows free movement of the operating table
- Unique shape designed to avoid pressure in the armpit, neck and genital area
- Regulatory compliance: CE, FDA, TGA, UKCA
- MD Class: I



## Technical Specifications

Chest	
Length	270 mm
Width	500 mm
Height	170 mm
Pelvis	
Length	356 mm
Width	500 mm
Height	170 mm
Weight capacity	160 kg
Product Code	CSM-2918-G
Manufacturer	Care Surgical Limited

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)