

Horizontal supporting roll 10369

Comfortable supporting roll for knee operations. Total length 454 mm, vertically adjustable 350 mm. Horizontally adjustable 175 mm. Turns horizontally 360 degrees. Pad size: length 300 mm, diameter Ø 133 mm.

Order no.

110369	Supporting Roll, horizontal 10369 incl. 10301 clamp
100 010 369	Supporting Roll, horizontal
100 010 301	Clamp 10301



Technical Specifications

Max. patient weight	250 kg
Safe working load (SWL)	40 kg
Weight	3 kg
Product Code	110369, 100010369

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com