

# K-Tek2603 Nissen Straps

K-Tek2603 Nissen Straps is used to prevent patient slipping during Steep Reverse Trendelenburg positioning. Mounts to any standard OR Rail. KTEK4700 Simple Blade Rail Clamps are required.



## Technical Specifications

Max. patient weight	227 kg
Product Code	100012822

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)