

# Lateral Arm Support, Multi Adjustable, SS4010

- Multi-adjustable arm support with ball joint
- Antimicrobial tightening lever
- Vertical tube length: 500 mm
- Viscoelastic cushion, 450 × 135 mm, with Velcro fastening
- Safe Working Load (SWL): 15 kg
- Recommended rail clamp: CL2110 (EU), CL2111 (US), CL2112 (UK)
- CE-marked medical device



## TECHNICAL SPECIFICATIONS

Product code	SS4010

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)