

## Lateral support, fixed SS1010

- fixed arthroscopy lateral support
- padding (50 mm) designed to prevent pressure and nerve injuries
- height adjustment range 150 mm
- pad size 100x200 mm
- Recommended clamps: CL2010/ CL2110



### TECHNICAL SPECIFICATIONS

Product code	SS1010

## Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)