


Shoulder support for lateral and side-lying position 1250P

- Pressure-relieving shoulder support cushion designed for surgical patients in lateral and side-lying positions
- Used in the operating room during hip and kidney procedures, as well as gynecological procedures
- The orthopedic shape provides a comfortable side-lying position with pressure relief
- Made of 100% SAFE ViscoPore® foam, offering high load-bearing capacity and long durability
- COMFOR™ OP covers feature a concealed plastic zipper and glued seams
- Strong, elastic, breathable, and antimicrobial.
- CE-marked Medical Device in accordance with MDR (EU) 2017/745.
- 85 x 50 x 13 cm (Length x Width x Height)
- Weight: 1,6 kg



Technical specifications

Manufacturer	TT SAFE MED ApS

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com

