


## Shoulder support, pair SS4040

- shoulder support keeps patient safely on operating table in Trendelenburg position
- soft, comfortable cushions prevent pressure and nerve injuries
- shoulder support attaches to the accessory rails
- Recommended clamp: CL2110



### TECHNICAL SPECIFICATIONS

Product code	SS4040

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)