

## CF-backrest 60502

Merivaara's Carbon fibre sections give surgeons full 360 deg. C-arm imaging during special procedures such as cardiovascular, orthopaedic and trauma operations. Attached to Promerix operating table's seat section normal or reverse mode.

Integrated carbon fibre EU-rails. Suitable clamps: CL2010, CL2020, CL2110, CL2120, CL2210.

Includes 80 mm VEF (Visco elastic foam) mattress, which is antistatic. Seams are welded.



## Technical Specifications

Length	1200 mm
Width (table top)	540 mm
Max. patient weight	275 kg
Product Code	100060502

### Contact us

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)