

CF-backrest 60502

Merivaara's Carbon fibre sections give surgeons full 360 deg. C-arm imaging during special procedures such as cardiovascular, orthopaedic and trauma operations. Attached to Promerix operating table's seat section normal or reverse mode.

Integrated carbon fibre EU-rails. Suitable clamps: CL2010, CL2020, CL2110, CL2120, CL2210.

Includes 80 mm VEF (Visco elastic foam) mattress, which is antistatic. Seams are welded.



Technical Specifications

Length	1200 mm
Width (table top)	540 mm
Max. patient weight	275 kg
Product Code	100060502

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com