

## Light non-motorized leg section 60203

Light leg section improves work ergonomics of the operating room staff. It is very light and easy to handle. Light weight leg section attaches to the operating table with Fluent locking mechanism. Light leg section for non-motorized legs, 100060203, is suitable for Smarter Practico operating table with non-motorized leg joints, Practico electro-mechanical and Practico Manual operating tables. To be used with Promerix operating table, 25- or 40- section is needed first to be attached to the seat section. Includes 80 mm thick, antistatic VEF (viscoelastic foam) mattress.



### Technical Specifications

Length	630-750 mm
Width (tabletop)	540 mm
Height	76 mm
Safe working load (SWL)	40 kg
Weight	4.9 kg
Product code	100060203

## **Contact us**

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)