

Leg rest with padding 19120

Radiotranslucent leg rest, curved support (carbon fibre) with support post and clamp, attaches to the extension bar. Rod diameter Ø 19 mm. Suitable with orthopedic extension device 19300.



Technical Specifications

Length	330 mm
Height	470 mm
Max. patient weight	225 kg
Safe working load (SWL)	29 kg
Weight	1.9 kg
Product Code	100019120

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com