

Support post for plaster cast 19132

For orthopedic extension device 19300. Rotates 360 degrees. Height adjustment 85–720 mm.

Rod diameter Ø 19 mm. Clamp 19134 to be ordered separately.



Technical Specifications

Length	780 mm
Height	400 mm
Max. patient weight	225 kg
Safe working load (SWL)	29 kg
Weight	3.9 kg
Product Code	100019132

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com