

Support post for plaster cast 19132

For orthopedic extension device 19300. Rotates 360 degrees. Height adjustment 85–720 mm.

Rod diameter \varnothing 19 mm. Clamp 19134 to be ordered separately.



Technical Specifications

| | |
|-------------------------|-----------|
| Length | 780 mm |
| Height | 400 mm |
| Max. patient weight | 225 kg |
| Safe working load (SWL) | 29 kg |
| Weight | 3.9 kg |
| Product Code | 100019132 |

Contact us

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com