

## Disposable sterile light handle cover

Disposable sterile light handle cover for Q-Flow surgical light and examination light Merilux X1

- Sterile lamp handles protection for single use
- Protection size fits various types of handles
- Lamp handle with protection 150 mm allows access the light switch on the handle, if applicable
- The perforated plastic collar acts both as plain protection against light and while also flexible holder
- The installed plastic cap creates a barrier between the operator and handle
- Sterile packed in peel pack, Ethylene oxide sterilization

A43460601 Set of 200 units





## Technical Specifications

Product Code	A43460600

## Связаться с нами

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)