


## Counter traction post 11438

Counter traction post 11438 for humerus, padded

- Requires a clamp 10300
- Pad size: length 300 mm, diam. 80 mm



### Technical Specifications

Width	365 mm
Height	639 mm
Weight	2 kg
Product code	110438, 100011438

### Kontakta oss

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)

