

# K-Tek Arm Boards

Integrated clamp snap locks onto rail at any position. Designed to prevent fluid ingress. Radiolucent, ultradurable strong material. Eliminates gear grinding and pinch hazards. One-hand adjustable. Order no. 100012818 K-TEK2005 Arm Board, EU-rail 100012819 K-TEK2000 Arm Board, US-rail 100012820 K-TEK2006 Arm Board, UK-rail



## Technical Specifications

Max. patient weight	227 kg
Weight	1.8 kg
Product code	100012818, 100012819, 100012820

## Kontakta oss

**Merivaara Corp.**  
Tarmontie 2-4  
15860 HOLLOLA  
FINLAND

+358 3 3394 611

[merivaara@merivaara.com](mailto:merivaara@merivaara.com)