

Tibia nailing kit for Flexion

Tibia Nailing Kit is meant for tibia nailing. It can be used with Flexion orthopedic extension device. Recommended to be used with Flexion trolley 519550, to be ordered separately.

Includes:

- Rail extension 19128
- Radial clamp 10300
- Nailing support A43592100
- Knee Crutch A43543900
- Traction device A43376200
- Stirrup Clamp 19114



Technical Specifications

Безопасная рабочая нагрузка	135 kg
код продукта	100519560

Kontakta oss

Merivaara Corp.
Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com