

SOLO

SOLO suspension solution for surgical lights and monitors

Modern operating rooms today are equipped with more and more equipment designed to provide surgical staff both the clinical images as well as critical video sources during the procedure. For the staff, functionality and ergonomics of suspension solutions that carry these vital monitors, is of very high importance.

Suspension solutions of Q-Flow product family has gone through a complete update, now being able to offer many new solutions that will enhance the ergonomics of any operating room.

[Read how Merivaara's R&D was challenged to think outside the box and what then happened.](#)

EXAMPLES OF AVAILABLE ARM OPTIONS:

- **SOLO Mobile**
- **SOLO Wall**
- **SOLO Ceiling** 50 mm, 65 mm or 110 mm ceiling tube



Technical Specifications

Product Code	520210

Kontakta oss

Merivaara Corp.

Tarmontie 2-4
15860 HOLLOLA
FINLAND

+358 3 3394 611

merivaara@merivaara.com